


Product Datasheet

Bcr (Phospho-Tyr360) antibody (orb106044)

Description	Rabbit polyclonal to Bcr (Phospho-Tyr360). Bcr is a serine/threonine and GTPase-activating protein (...)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200, IF/ICC:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of COS7 whole cell...