


Product Datasheet

ATF2 (Phospho-Ser112 or 94) antibody (orb106060)

Description	Rabbit polyclonal to ATF2 (Phospho-Ser112 or 94). This gene encodes a transcription factor that is a...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
MDA-MB-
435 whole...