


Product Datasheet

ATF2 (Phospho-Thr73 or 55) antibody (orb106063)

Description	Rabbit polyclonal to ATF2 (Phospho-Thr73 or 55). This gene encodes a transcription factor that is a ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 whole cell ...