


Product Datasheet

GATA1 (Phospho-Ser310) antibody (orb106066)

Description	Rabbit polyclonal to GATA1 (Phospho-Ser310). GATA1 Transcriptional activator which probably serves a...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of COS7 whole cell...