



Product Datasheet

ERK1 (phospho-Thr202/Tyr204) antibody (orb10607)

Description	ERK1 (phospho-Thr202/Tyr204) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from rat ERK1 around the phosphorylation site of Thr201/204 (FL(p-T)E(p-Y)VA)
Target	phospho-Erk1 (Thr202 + Tyr204)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	Q63844
Dilution Range	IHC-P=1:100-500, Flow-Cyt=2ug/Test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brain...



Flow cytometric analysis of U251 cell us...