


Product Datasheet

PECAM-1 (Phospho-Tyr713) antibody (orb106070)

Description	Rabbit polyclonal to PECAM-1 (Phospho-Tyr713). PECAM-1 is a cell adhesion molecule expressed on plat...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IF/ICC:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat whole ce...