

Product Datasheet

Stathmin 1 (Phospho-Ser24) antibody (orb106073)

Description Rabbit polyclonal to Stathmin 1 (Phospho-Ser24).
This gene belongs to the stathmin family of genes.
...

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated


Tested Applications IHC-P, WB

Note For research use only

Clonality Polyclonal

Dilution Range WB:1:500-1:2000, IHC-P:1:50-1:200

Expiration Date 12 months from date of receipt.



Western
blot
analysis of
Jurkat
whole ce...