

Product Datasheet

Stathmin 1 (Phospho-Ser37) antibody (orb106074)

Description Rabbit polyclonal to Stathmin 1 (Phospho-Ser37). This gene belongs to the stathmin family of genes. ...

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated


Tested Applications IHC-P, WB

Note For research use only

Clonality Polyclonal

Dilution Range WB:1:500-1:2000, IHC-P:1:50-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of Jurkat whole ce...