


Product Datasheet

Connexin 43 (Phospho-Ser367) antibody (orb106077)

Description	Rabbit polyclonal to Connexin 43 (Phospho-Ser367). Gap junction protein, alpha 1 is a member of the ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
K562 whole
cell...