

Product Datasheet

GABA-RB (Phospho-Ser434) antibody (orb106083)

Description Rabbit polyclonal to GABA-RB (Phospho-Ser434). GABRB1 GABA, the major inhibitory neurotransmitter in...

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated


Tested Applications IHC-P, WB

Note For research use only

Clonality Polyclonal

Dilution Range WB:1:500-1:2000, IHC-P:1:50-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of COS7 whole cell...