


Product Datasheet

NF-kappaB p105/p50 (Phospho-Ser337) antibody (orb106092)

Description	Rabbit polyclonal to NF-kappaB p105/p50 (Phospho-Ser337). NFkB-p105 a transcription factor of the nu...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of MDA-MB-435 whole cell lysate.