


Product Datasheet

p70 S6 kinase (Phospho-Ser424) antibody (orb106097)

Description	Rabbit polyclonal to p70 S6 kinase (Phospho-Ser424). This gene encodes a member of the RSK (ribosoma...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 whole cell ...