


Product Datasheet

Cofilin (Phospho-Ser3) antibody (orb106099)

Description	Rabbit polyclonal to Cofilin (Phospho-Ser3). cofilin 1 a cytoskeletal protein that controls actin de...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of COLO205 whole c...