


Product Datasheet

CDC25A (Phospho-Ser75) antibody (orb106113)

Description	Rabbit polyclonal to CDC25A (Phospho-Ser75). CDC25A is a member of the CDC25 family of phosphatases....
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of A2780 whole cel...