

Product Datasheet

CDC25A (Phospho-Ser124) antibody (orb106115)

Description Rabbit polyclonal to CDC25A (Phospho-Ser124). CDC25A is a member of the CDC25 family of phosphatases...

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated


Tested Applications WB

Note For research use only

Clonality Polyclonal

Dilution Range WB:1:500-1:2000

Expiration Date 12 months from date of receipt.



Western blot analysis of 293 whole cell ...