

Product Datasheet

Catenin-beta (Phospho-Ser33) antibody (orb106126)

Description	Rabbit polyclonal to Catenin-beta (Phospho-Ser33). Beta-catenin is an adherens junction protein. Adh...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.