


Product Datasheet

iNOS (Phospho-Tyr151) antibody (orb106129)

Description	Rabbit polyclonal antibody to iNOS (Phospho-Tyr151)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
whole cell
lysa...