



Product Datasheet

IRS1 (phospho-Ser312) antibody (orb106132)

Description	IRS1 (phospho-Ser312) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human IRS1 around the phosphorylation site of Ser312 (AT(p-S)PA)
Target	phospho-IRS1 (Ser312)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	132 kDa
Uniprot ID	P35568
Dilution Range	IHC-P=1:100-500, Flow-Cyt=0.2µg/Test
Expiration Date	12 months from date of receipt.



Flow cytometry of HeLa(blue line).Primar...



IHC-P analysis of rat pancreas using pho...