

Product Datasheet

IRS-1 (Phospho-Ser639) antibody (orb106134)

Description	Rabbit polyclonal to IRS-1 (Phospho-Ser639). IRS-1 is an adaptor protein that is one of the major su...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.