


Product Datasheet

VEGFR2 (Phospho-Tyr951) antibody (orb106137)

Description	Rabbit polyclonal antibody to VEGFR2 (Phospho-Tyr951)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of SK-OV3 whole ce...