


Product Datasheet

MEK2 (Phospho-Thr394) antibody (orb106138)

Description	Rabbit polyclonal to MEK2 (Phospho-Thr394). MEK2 a dual-specificity protein kinase of the STE7 kinas...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of ovarycancer who...