


## Product Datasheet

### STAT3 (Phospho-Ser727) antibody (orb106146)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to STAT3 (Phospho-Ser727). The protein encoded by this gene is a member of the STA... |
| <b>Species/Host</b>        | Rabbit                                                                                                  |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                            |
| <b>Tested Applications</b> | IHC-P, IP, WB                                                                                           |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                              |
| <b>Dilution Range</b>      | WB:1:500-1:2000, IHC-P:1:50-1:200                                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                         |



Western blot analysis of HeLa whole cell...