


Product Datasheet

STAT1 (Phospho-Tyr701) antibody (orb106150)

Description	Rabbit polyclonal to STAT1 (Phospho-Tyr701). The protein encoded by this gene is a member of the STA...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF7 whole cell...