

Product Datasheet

SEK1/Mkk4 (Phospho-Thr261) antibody (orb106165)

Description Rabbit polyclonal to SEK1/Mkk4 (Phospho-Thr261). MKK4 dual specificity kinase of the STE7 family tha...

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated


Tested Applications IHC-P, WB

Note For research use only

Clonality Polyclonal

Dilution Range WB:1:500-1:2000, IHC-P:1:50-1:200

Expiration Date 12 months from date of receipt.



Western
blot
analysis of
293 whole
cell ...