

Product Datasheet

p27 kip1 (Phospho-Thr187) antibody (orb106166)

Description Rabbit polyclonal to p27 kip1 (Phospho-Thr187). This gene encodes a cyclin-dependent kinase inhibito...

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated


Tested Applications WB

Note For research use only

Clonality Polyclonal

Dilution Range WB:1:500-1:2000

Expiration Date 12 months from date of receipt.



Western blot analysis of HeLa whole cell...