


Product Datasheet

p27 kip1 (Phospho-Ser10) antibody (orb106167)

Description	Rabbit polyclonal to p27 kip1 (Phospho-Ser10). This gene encodes a cyclin-dependent kinase inhibitor...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of A2780 whole cel...