


Product Datasheet

GSk3 alpha (Phospho-Ser21) antibody (orb106172)

Description	Rabbit polyclonal to GSk3 alpha (Phospho-Ser21). GSK3A a proline-directed protein kinase of the GSK ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
ovarycancer
who...