


Product Datasheet

CCR5 (Phospho-Ser336) antibody (orb106175)

Description	Rabbit polyclonal to CCR5 (Phospho-Ser336). CCR5 a 7-transmembrane G-linked receptor for a number of...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat whole ce...