



Product Datasheet

EpCAM antibody (orb10618)


Description	EpCAM antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human EpCAM. Please contact us for the exact immunogen sequence. The peptide is available as orb374722.
Target	EpCAM
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	35 kDa
Uniprot ID	Q3T0L5
NCBI	002354 , 21 , 002345
Entrez	4072
Dilution Range	IHC-P: 1:200
Expiration Date	12 months from date of receipt.



IHC-P staining of human endometrial cancer tissue.



IHC-P image of mouse lung tissue using a rabbit anti-human EpCAM antibody.



Immunohistochemical staining of paraffin-embedded tissue.