


Product Datasheet

CDk5 (Phospho-Tyr15) antibody (orb106197)

Description	Rabbit polyclonal to CDk5 (Phospho-Tyr15). CDK5 is a protein kinase of the CDK family. Unlike other ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of COS7 whole cell...