


Product Datasheet

MEF2A (Phospho-Thr312) antibody (orb106205)

Description	Rabbit polyclonal to MEF2A (Phospho-Thr312). MEF2A a myocyte-specific enhancing transcription factor...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH-3T3 whole c...