



Product Datasheet

NFkB p65 (phospho-Ser276) antibody (orb106211)


Description	NFkB p65 (phospho-Ser276) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human NFkBp65 around the phosphorylation site of Ser276 (RP(p-S)DR)
Target	Phospho-NFkB p65 (Ser276)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61 kDa
Uniprot ID	Q04206
Dilution Range	IHC-P=1:100-500, ICC=1:100, Flow-Cyt=0.2µg /Test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse lymph node tissue.



Immunohistochemical staining of mouse lymph node tissue.



Immunohistochemical staining of rat spleen tissue.