


Product Datasheet

NFkB p65 (Phospho-Thr435) antibody (orb106216)

Description	Rabbit polyclonal antibody to NFkB p65 (Phospho-Thr435)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of COS7 whole cell...