


Product Datasheet

Dynamin-1 (Phospho-Ser774) antibody (orb106220)

Description	Rabbit polyclonal to Dynamin-1 (Phospho-Ser774). Dynamins represent one of the subfamilies of GTP-bi...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
Mouse
brain tis...