


Product Datasheet

PkC zeta (Phospho-Thr560) antibody (orb106226)

Description	Rabbit polyclonal to PkC zeta (Phospho-Thr560). Protein kinase C (PKC) zeta is a member of the PKC f...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of COS7 whole cell...