


Product Datasheet

4E-BP1 (Phospho-Thr45) antibody (orb106242)

Description	Rabbit polyclonal to 4E-BP1 (Phospho-Thr45). 4E-BP1 binds to eIF4E, preventing its assembly into the...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of MDA-MB-435 whole cell lysate.