


Product Datasheet

PkD1/PkCmicro (Phospho-Tyr463) antibody (orb106252)

Description	Rabbit polyclonal to PkD1/PkCmicro (Phospho-Tyr463). Members of the protein kinase C (PKC) family fu...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
HepG2
whole cel...