


Product Datasheet

p38 MAPK (Phospho-Thr180) antibody (orb106259)

Description	Rabbit polyclonal to p38 MAPK (Phospho-Thr180). The protein encoded by this gene is a member of the ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
HeLa whole
cell...