


Product Datasheet

DRP-2 (Phospho-Thr514) antibody (orb106260)

Description	Rabbit polyclonal to DRP-2 (Phospho-Thr514). CRMP-2 is an enzyme with dihydropyrimidinase activity. ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200, IF/ICC:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH-3T3 whole c...