


Product Datasheet

PDCD4 (Phospho-Ser67) antibody (orb106264)

Description	Rabbit polyclonal to PDCD4 (Phospho-Ser67). This gene encodes a protein localized to the nucleus in ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 whole cel...