


Product Datasheet

FRS2 (Phospho-Tyr436) antibody (orb106266)

Description	Rabbit polyclonal to FRS2 (Phospho-Tyr436). FRS2 is an adaptor protein involved in fibroblast growth...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of COS7 whole cell...