


Product Datasheet

BLNk (Phospho-Tyr96) antibody (orb106267)

Description	Rabbit polyclonal to BLNk (Phospho-Tyr96). BLNK an adaptor protein that bridges the B-cell receptor...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of COLO205 whole c...