



Product Datasheet

GTPase Activating Protein (Phospho-Ser387) antibody (orb106281)


Description	Rabbit polyclonal to GTPase Activating Protein (Phospho-Ser387). Rho GTPases control a variety of ce...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC 1:100-1:500, ELISA(peptide): 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of Heat-shock trea...



Western blot analysis of COS7 cells lysa...



Immunohistochemical staining human breas...