


Product Datasheet

HDAC5 (Phospho-Ser498) antibody (orb106286)

Description	Rabbit polyclonal to HDAC5 (Phospho-Ser498). HDAC9 a transcriptional regulator of the histone deacet...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH-3T3 whole c...