

Product Datasheet

Ikk-gamma (Phospho-Ser31) antibody (orb106289)

Description Rabbit polyclonal to Ikk-gamma (Phospho-Ser31). Familial incontinentia pigmenti (IP) is a genodermat...

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated


Tested Applications IHC-P, WB

Note For research use only

Clonality Polyclonal

Dilution Range WB:1:500-1:2000, IHC-P:1:50-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of 293 whole cell ...