


Product Datasheet

SARS-CoV-2 (2019-nCoV) Nucleocapsid antibody (orb1063083)

Description	SARS-CoV-2 (2019-nCoV) Nucleocapsid antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protein (full length)
Target	SARS-CoV-2 (2019-nCoV) Nucleocapsid
Preservatives	1*PBS, 0.03% Proclin300.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	46 kDa
Uniprot ID	P0DTC2
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis using Anti-SARS-CoV-2 Nucleocapsid antibody (orb1063083).