


Product Datasheet

Dab1 antibody (orb106310)

Description	Rabbit polyclonal to Dab1. DAB1 an adaptor molecule functioning in neural development. The laminar o...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000, IF/ICC:1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 whole cel...