



Product Datasheet

Factor VIII antibody (orb10633)


Description	Factor VIII antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from rat factor VIII (401-500/2258 aa)
Target	factor VIII
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	267 kDa
Uniprot ID	P00451
Dilution Range	IHC-P=1:100-500, ICC=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IF staining mouse c57 cells (OPG antibod...



Immunohistochemical staining of paraffin...



IHC-p of pig tissue (orb10633 at 1:200)