

Product Datasheet

SARS-CoV-2 (2019-nCoV) Nucleocapsid, WT & Ormicron Recognized antibody (orb1063822)

Description SARS-CoV-2 (2019-nCoV) Nucleocapsid, WT & Ormicron Recognized antibody

Species/Host Mouse

Reactivity Virus

Conjugation Unconjugated

Tested Applications WB

Immunogen Recombinant SARS-CoV-2 Nucleocapsid Protein (1-419/419 aa)

Target SARS-CoV-2 (2019-nCoV) Nucleocapsid, WT & Ormicron Recognized

Preservatives 0.01M PBS(pH7.4)

Form/Appearance Liquid

Storage Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note For research use only

Isotype IgG1


Clonality Monoclonal

MW 46 kDa


Uniprot ID **P0DTC2**

Dilution Range WB=1:500-2000

Expiration Date 12 months from date of receipt.



Western blot analysis using Anti-SARS-CoV-2



Western blot analysis using Anti-SARS-CoV-2