

Product Datasheet

SARS-CoV-2 (2019-nCoV) Spike RBD (WT & Omicron) antibody (orb1063825)

Description	SARS-CoV-2 (2019-nCoV) Spike RBD (WT & Omicron) antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant SARS-CoV-2 Spike S1 Protein (14-685/1273 aa)
Target	SARS-CoV-2 (2019-nCoV) Spike RBD (WT & Omicron)
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Concentration	Lot Dependent
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	140 kDa
Uniprot ID	P0DTC2
Expiration Date	12 months from date of receipt.